

Ureteral Dilation Balloon Catheters

M-range for rigid ureterscope

soft tip guidewire mounted on the dilation balloon

balloon go through working channel of rigid endoscope, under direct vision of endoscope when placing dilation balloon

hydrophilic Coating

N-range for guidewire (standard)

RO marker

precise positioning scale catheter

order information

Product Code	Description	Balloon diameter x length	Catheter diameter x length
wireguide type(N-range)			
IVX-UB-UY07060	dilation balloon catheter, include inflator	7mm x 60mm	6Fr x 700mm
ureterscope type(M-range)			
IVX-UB-MY04040	dilation balloon catheter,	4 mm x 40 mm	2.4F x 1500mm
IVX-UB-MY05060	soft-tip guidewire mounted on balloon,	5 mm x 60 mm	2.4F x 1500mm
IVX-UB-MY07060	include inflator	7 mm x 60 mm	2.4F x 1500mm