


# Single Use Grasping Forceps




## order information

| Product Code  | Description                              | Specification |
|---------------|------------------------------------------|---------------|
| IVX-GF1-1808A | Sheath OD:1.9mm,<br>working length: 80cm | 1.9mmx 80cm   |

## INNOVEX MEDICAL CO.,LTD.

Add: 2nd Floor, No.17 Building, 315 Qingda Rd, Pudong,  
201201 Shanghai, P.R. China  
Tel: +0086 021 50791830  
Web: <http://www.innovexmed.com.cn> Zip: 201201  
E-mail: [global@innovexmed.com.cn](mailto:global@innovexmed.com.cn)


# DISPOSABLES ENDOUROLOGY



[www.innovexmed.com.cn](http://www.innovexmed.com.cn)

# Ureteral Access Sheaths


negative pressure & tip flexible



## order information

| Product Code | Description        | Specification |
|--------------|--------------------|---------------|
| IVX-NS-1040  |                    | 10Fr x 40cm   |
| IVX-NS-1050  | tip flexible -     | 10Fr x 50cm   |
| IVX-NS-1140  | negative pressure  | 11Fr x 40cm   |
| IVX-NS-1150  | access sheath with | 11Fr x 50cm   |
| IVX-NS-1240  | scale, hydrophilic | 12Fr x 40cm   |
| IVX-NS-1250  | sheath and dilator | 12Fr x 50cm   |

# Ureteral Access Sheaths




## order information

| Product Code | Description        | Specification |
|--------------|--------------------|---------------|
| IVX-UA-1035  |                    | 10Fr x 35cm   |
| IVX-UA-1045  |                    | 10Fr x 45cm   |
| IVX-UA-1135  |                    | 11Fr x 35cm   |
| IVX-UA-1145  | access sheath with | 11Fr x 45cm   |
| IVX-UA-1235  | scale, hydrophilic | 12Fr x 35cm   |
| IVX-UA-1245  | sheath and dilator | 12Fr x 45cm   |

# Urological Guidewires

safe, smooth, easy navigation



USO - PU guidewire

BSO - hybrid guidewire (PTFE zebra)  
(Structure: floppy black PU head (Hydrophilic+ Tungsten provides radiopacity.) + PTFE Zebra, nitinol core)

TSO - ball tip PTFE guidewire

NSO - hybrid guidewire  
(stainless steel wire wounded sheath)

## order information

| Product Code   | Description                                                                                                  | Specification            |
|----------------|--------------------------------------------------------------------------------------------------------------|--------------------------|
| IVX-UG-3515USO | full length black PU (hydrophilic), nitinol core                                                             | .035in x 1500mm straight |
| IVX-UG-3215BSO | floppy black PU head (Hydrophilic+ Tungsten provides radiopacity.) + PTFE Zebra, nitinol core                | .032in x 1500mm straight |
| IVX-UG-3515BSO |                                                                                                              | .035in x 1500mm straight |
| IVX-UG-3215TSO | ball tip + PTFE Zebra, nitinol core                                                                          | .032in x 1500mm straight |
| IVX-UG-3515TSO |                                                                                                              | .035in x 1500mm straight |
| IVX-UG-3515NSO | floppy black PU head (hydrophilic) + stainless steel wire wounded sheath (PTFE powder coating), nitinol core | .035in x 1500mm straight |


# Urological Stone Extraction Baskets

4 wire tipless basket

3 wire basket

multi-layer sheath provides exceptional strength and flexibility

ergonomic handle



## order information

| Product Code | Description                                                   | Specification (opening x length) |
|--------------|---------------------------------------------------------------|----------------------------------|
| IVX-US-F1110 | nitinol stone retrieval basket, tipless, 1.9Fr sheath, 4 wire | 11mm x 1300mm                    |
| IVX-US-T0810 | nitinol stone retrieval basket, 1.9Fr sheath, 3 wire          | 8mm x 1300mm                     |
| IVX-US-T1110 |                                                               | 11mm x 1300mm                    |

## Ureteral Dilation Balloon Catheters

### M-range for rigid ureteroscope

soft tip guidewire mounted on the dilation balloon

balloon go through working channel of rigid endoscope, under direct vision of endoscope when placing dilation balloon

hydrophilic coating

### N-range for guidewire (standard)

RO marker

#### order information

| Product Code                       | Description                                    | Balloon diameter x length | Catheter diameter x length |
|------------------------------------|------------------------------------------------|---------------------------|----------------------------|
| <b>wireguide type (N-range)</b>    |                                                |                           |                            |
| IVX-UB-UY07060                     | dilation balloon catheter, include inflator    | 7mm x 60mm                | 6Fr x 700mm                |
| <b>ureteroscope type (M-range)</b> |                                                |                           |                            |
| IVX-UB-MY04040                     | dilation balloon catheter, soft-tip            | 4mm x 40mm                | 2.4Fr x 1500mm             |
| IVX-UB-MY05060                     | guidewire mounted on balloon, include inflator | 5mm x 60mm                | 2.4Fr x 1500mm             |
| IVX-UB-MY07060                     | balloon, include inflator                      | 7mm x 60mm                | 2.4Fr x 1500mm             |

precise positioning scale catheter

## Ureteral Stent

#### order information

| Product Code     | Description                                       | Specification |
|------------------|---------------------------------------------------|---------------|
| IVX-UT-01-472610 | single layer PU ureteral stent with push catheter | 4.7Fr x 26cm  |
| IVX-UT-01-602610 |                                                   | 6Fr x 26cm    |

## Urological Sampler Sets

sheath connection tube

ON/OFF valve

metal mesh filter

bottle with hook

less volume  
less weight

#### order information

| Product Code | Description                                 | Specification |
|--------------|---------------------------------------------|---------------|
| IVX-NB-150   | stone fragment collection with ON/OFF valve | 70mL          |