



EOS-H-01



EA100


**VIDEO
PROCESSOR**

6.3Fr

SINGLE USE FLEXIBLE
URETERORENOSCOPE

US23K-10-EU/US

**NEW
DESIGN**



9Fr ureteral access sheath

Slim without sacrificing performance
of **1.05mm** round working channel

**Shanghai AnQing
Medical Instrument Co., Ltd**

3&4 Floor, No.2 Building, 366 Huiqing Rd,
East Zhangjiang High-Tech Park,
201201 Shanghai, China
Tel: +86 21 61117375
Fax: +86 21 61117374
Email: global@innovexmed.com.cn

BRIGHT & CLEAR

HIGH IMAGE QUALITY

6.3Fr

SINGLE USE FLEXIBLE
URETERORENOSCOPE



Superior deflection performance

Excellent maneuverability & torque stability

Cable position at lower part of handle

Rotation function improve torque stability further

Nominal width of insertion section 6.3Fr

Insertion section length 670mm

Photo / video bottom

Lever lock

275°

Dual deflection



Inner layer PTFE improves smoothness for accessories



Field of view: 130°±15%
Depth of field: 3mm ~ 100mm



1.05mm round working channel



Super mini CMOS 160K



Dual LED illumination