



EA100

Video Processor



EOS-H-01

Shanghai AnQing Medical Instrument Co., Ltd


3&4 Floor, No.2 Building, 366 Huiqing Rd,
East Zhangjiang High-Tech Park,
201201 Shanghai, China

Tel: +86 21 61117375
Fax: +86 21 61117374
Email: global@innovexmed.com.cn

7.5FR

SINGLE USE FLEXIBLE URETERORENOSCOPE

• US27K-12-EU/US •



9.5Fr ureteral
access sheath

Slim without sacrificing
performance of **1.2mm**
round working channel



7.5FR

SINGLE USE FLEXIBLE URETERORENOSCOPE

Photo / video button

Lever lock

Cable position at
lower part of handle

Rotation function improve
torque stability further

Superior deflection
performance

Excellent maneuverability
& torque stability

Nominal width of
insertion section 7.5Fr

Insertion section
length 670mm

BRIGHT & CLEAR

HIGH IMAGE QUALITY

Order Information:

20 units / transportation carton = 2 units / carton x 10 cartons / transportation carton



275°
Dual deflection



1.2mm round working channel



Inner layer PTFE improves
smoothness for accessories



CMOS 160K



Dual LED illumination



Field of view: 110°±15%
Depth of field: 3mm ~ 100mm